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by destroying glass, cutting or tearing the primary, secondary and final packaging (box) and then they must be clearly and visibly identified, indicating your situation, that is, UNFIT FOR USE. The methods and procedures for disposing of the medical device must guarantee its complete lack of characterization, preventing any possibility of reuse. The mischaracterization of the medical product is the sole responsibility of the medical surgeon or dentist who is using the product, as well as the disposal methods and procedures used.

CARE AND PRECAUTIONS

For safety and effectiveness in the use of the GenMix graft, the following precaution and care measures must be observed:
- Only qualified professionals with knowledge of the surgical techniques and procedures necessary for proper choice and application should use the product.

- The medical surgeon or dentist must be familiar with the procedures and recommendations for implanting GenMix (see USE INFORMATION).

- Only specific sterile surgical instruments should be used.

- The surgeon must observe the conditions of preparation and asepsis of the area and instruments, and vestments of the team.
- As with all surgical procedures, care should be taken when treating medically compromised patients such as patients on steroid therapy for a long time, or patients on anticoagulant therapy.

GenMix is a SINGLE USE PRODUCT and must not be re-sterilized. Open products must be discarded.

WARNINGS

The surgeon must fully understand not only the medical and surgical aspects of the medical device.

The surgeon must take care to alert or test the patient for sensitivity to products based on bovine derivatives.

Smoking, drug users, alcoholics or malnourished patients should be warned about the increased incidence of bone non-union.

Note: All bovine bone substitutes from the Genius line are manufactured with animal bones from the Brazilian herd. According to the geographical risk framework issued by the International Zoonasitary Code and the Scientific Steering Comitee of the European Community (SSC EC August 2005), Brazil is free of Bovine Spongiform Encephalopathy (BSE). Although through decree 516/97, Brazil declares itself free from Bovine Enpongiform Encephalopathy, although the processing to which the products are submitted is known to be effective in inactivating the causative agent of the EEBs, although the animals used for the production of the Genius line are registered in the Brazilian system of bovine and buffalo traceability - SISBOV, like all products of bovine origin, even if remote, there is a risk of transmission of BSE. In cases of appearance of symptoms of nvCDJ the health professional must notify the health authority (tecnovigilância. @Anvisa.gov.br)

POSSIBLE ADVERSE EFFECTS

No adverse reactions have been reported for the GenMix product.

Possible complications such as swelling, infections, pain and complications associated with anesthesia may occur because they are common to implant surgery. As with any type of implant surgery, the patient may experience discomfort in the first postoperative days.

In addition to the aforementioned, normal reactions in periodontal surgeries may be found such as: swelling of intraoral tissue, thermal sensitivity, gum recession, some loss of height of crystal bone, etc.

In the clinical studies evaluating the GenMix product, no immune reactions or hypersensitivity reactions were found.

The surgeon must be fully informed about the medical and surgical aspects of the medical device.

The surgeon must take care to alert or test the patient for sensitivity to products based on bovine derivatives.

Smoking, drug users, alcoholics or malnourished patients should be warned about the increased incidence of bone non-union.

SPECIAL CARE AND CLARIFICATION ON MEDICAL PRODUCT STORAGE AND TRANSPORTATION

The packaging must be intact upon receipt.

Open packaging should not be used (do not use the medical product if the packaging is broken);

The product must be transported and stored at room temperature. Avoid humidity and excessive heat;

The storage location of the medical device must be clean, dry and bright in order to maintain the ideal conditions for storage and transport, as well as its physical and chemical integrity;

The carrier must be informed of the content and delivery time. The precautions to be taken during transport are demonstrated on the packaging, thus ensuring protection of the medical product from dispatch to delivery to the customer.

CONTRAINDICATIONS

The GenMix product is not designed for any use other than that indicated. Do not use the product in the presence of any of the contraindications listed below.

The use of GenMix is not recommended in patients with the following clinical conditions:

- Uncontrolled metabolism disorders
- Hypersensitivity or allergic reactions to collagen
- Hypersensitivity or allergic reactions to bovine derivatives
- Generalized renal dysfunction
- Severe liver disease
- Bone immaturity;
- Osteopenia or severe osteoporosis;
- Presence or recent history of systemic or localized infection (contaminated wounds in or near the implant)
- Alcoholism or drug abuse;
- Pregnancy and lactation (at the doctor’s discretion);
- Fever or leukocytosis;

The use of GenMix is not recommended in situations where periodontal surgery is not indicated

GenMix should not be used in the presence of tumors or a history of isolated malignancy or associated with deproteinized or demineralized bovine bone matrix or synthetic or auto or allogeneic grafts.

The use of GenMix is not recommended in any medical or surgical condition that may compromise the success of the surgical procedure;

STERILITY

The GenMix medical product is marketed sterile.

Leftover product should not be reused.

GenMix cannot be re-sterilized.

WARRANTY

The guarantee will be applied to medical products as long as it is used in accordance with the previous sections of this instruction for use: SUPPORT MATERIALS, INDICATIONS, USE INFORMATION, CONTRAINDICATIONS, CAUTIONS AND PRECAUTIONS, WARNING, SPECIAL CARES AND CLARIFICATIONS ABOUT THE STORAGE AND TRANSPORT OF MEDICAL PRODUCT and STERILITY.

COMPLAINT

Any customer or user of these medical products who wishes to claim, if he is not satisfied with the services and / or products offered by BAUMER, such as: Identification, traceability, security, effectiveness, performance, should contact the authorized BAUMER distributor.

Any problems or suspected problems with related medical products should be reported immediately to BAUMER or your authorized BAUMER distributor.

Any serious adverse effects that affect the health or safety of the patient or user, such as malfunction or damage to medical products, serious problems or death caused by, or associated with, the use of these products, must be reported to BAUMER.

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Baumer S.A