# Genius

Pratix

## Membrana Regeneradora

### Membrana Biológica Sintética - Absorvível Pratix

#### DESCRIÇÃO

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### Membrana Biológica Sintética - Absorvível Pratix

#### DESCRIÇÃO

Membrana Biológica Sintética – Absorvível (Pratix) é a synthetic, implantable, resorbable membrane designed for use as a biological barrier by dental professionals in guided bone or tissue regeneration procedures. It acts as a mechanical barrier that prevents the migration of epithelial tissue and gingival connective tissue to the defect area, maintaining space between the membrane and the bone, allowing cell repopulation, bone maturation and defect regeneration.

Membrana Biológica Sintética – Absorvível (Pratix) is a biological barrier composed of poly (lactic acid-co-glycolic acid) PLGA copolymers. These copolymers are degraded in vivo by hydrolysis when in contact with body fluids generating products that will be metabolized by the body resulting in CO2 and H2O. As it is resorbable, the use of Membrana Biológica Sintética – Absorvível (Pratix) eliminates the need for a second surgery to remove the membrane.

Membrana Biológica Sintética – Absorvível (Pratix) is an easy to handle and highly flexible membrane, with a controlled degradability rate in a physiological environment. The degradation time and consequent resorption can vary between 8 to 10 weeks.

Membrana Biológica Sintética – Absorvível (Pratix) is presented as a rectangular membrane and disk in different sizes for better adaptation to the defects to be corrected.

Membrana Biológica Sintética – Absorvível (Pratix) is a STERILE product, sterilized with GAMMA RADIATION, SINGLE USE PRODUCT, therefore it is PROHIBITED TO REPROCESS.


Figura 1 - Macroscopic aspect of the Synthetic Biological Membrane - Absorbible (Pratix), in a rectangular and disc shape.

#### FUNDAMENTALS OF PRODUCT ACTION

One of the goals of periodontal therapy is to correct defects or anatomical problems caused by the disease process. Various surgical techniques have been suggested in an attempt to regenerate, partially or totally, the tissues lost in advanced periodontal disease.

In particular, techniques and tools have evolved to inhibit the growth of adjacent tissues that can replace lost or harmful tissue. In recent years, the principle of guided tissue regeneration (RTG) has been developed and improved for the treatment and regeneration of periodontal tissues lost due to inflammatory processes in the periodontium and the regeneration of bone tissue in edentulous areas that would be subjected to the placement of osteo-integrable metal implants. Conceptually, regeneration is a situation where the repair process, after tissue injury, culminates in the complete restoration of the morphological and functional aspects of the tissue (TEN CATE, 1994 apud Taga, 2004).

The principle of RTG is to use a physical barrier to create a suitable environment, protected from the invasion of competing, less differentiated and highly proliferative cells from neighboring tissues, mainly fibroblasts and epithelial cells, and to promote the proliferation of cells originally characteristic of the region for regeneration. of the natural tissues of the area (NYMAN et al., 1980 apud Taga, 2004).

From the 1980s, biocompatible, absorbable and non-absorbable membranes were proposed and commercialized based on the theory that cells from the periodontal ligament are the only ones with the potential for differentiation in cement cells, with the membranes acting as a barrier preventing the migration of gingival fibroblasts

and epithelial cells over the naked root, allowing selective relocking by the cells of the remaining periodontal ligament. Currently, the RTG technique has also been employed in the treatment of angular bone defects around osseous integrated implants.

Two classes of barriers are used in RTG procedures: absorbable and non-absorbable. The indication of which type of membrane should be used, is directly linked to the type of defect chosen, the cost-benefit ratio and the patient's availability or interest in performing reentry surgery to remove the non-absorbable membrane. According to the proposal by Wang et al 1998 apud Taga 2004, such barriers must have the following characteristics: a) biocompatibility b) permeability, allowing the diffusion of plasma and nutrients, but not the passage of cells c) act as a physical support to the adjacent soft tissue, preventing its collapse in the space with the blood clot, where the cells of the original tissue of the wound may migrate to, d) provide protection of the simple vascular network while the clot undergoes an organization process e) perform exclusion of competing cells f) prevent the formation of a tissue scar and g) to be easy to handle. The advantages of an absorbable material over non-absorbable material would be:

- 1 - The elimination of the second surgery to remove the barrier, which allows for a reduction in operating time and, consequently, costs; b) reduction of the patient's morbidity; c) increase in the general acceptance of the patient for RTG and Osteopromotion procedures; and d) reduction of the risks of insertion loss obtained by regeneration due to reentry surgery.
- 2 - Biologically absorbable materials have the potential to: a) resist or prevent bacterial colonization; b) increase tissue coverage or protection and reduce the exposure of the barrier; and c) have a greater integration with the host tissue.

Pratix is a resorbable membrane composed of copolymers of polylactic acid and polylglycolic acid. The production method of Pratix allows to obtain a biocompatible product and with absorption properties by the tissues of the host ranging around 8 and 10 weeks (depending on the thickness), as has been demonstrated in animal studies and when implanted in humans.

The membranes of these copolymers are absorbed by hydrolysis and once re-absorbed, they are replaced by endogenous scar tissue and not any type of tissue or material foreign to the organism has been detected.

Its final shape has ideal porosity that allows a selective action of filtering cells and nutrients to the implant site where bone remodeling will occur. The flexibility of these membranes is of fundamental importance in adapting them to the anatomical structures where they will be inserted during the operation. Pratix acts as a mechanical barrier in the Guided Tissue Regeneration (RTG) technique, which is fundamental in the use of a passive barrier to prevent the migration of epithelial tissue and gingival connective tissue to the bone defect during remodeling.

When these barriers are placed between the flap and the root surface, contact of the connective and epithelial tissue with that surface is avoided. This allows cells derived from the periodontal ligament to be selectively induced on the root surface, favoring regeneration than periodontal tissue. Once the Membrana Biológica Sintética – Absorvível (Pratix) is placed, it will act as a biological barrier preventing the invagination of soft tissue into the implant area and interfering with bone regeneration and repair.

The absorption speed of the membranes can be influenced by local characteristics, such as vascularization of the region, thickness of the mucoza, degree of inflammation resulting from the surgery itself, changes in the pH of the tissue involved in the repair, local infections and partial exposures to the oral environment.

#### INDICAÇÃO

Membrana Biológica Sintética – Absorvível (Pratix) is an absorbable, implantable polylactic glycolic membrane indicated for use in tissue guided regeneration in the covering of other materials in periodontal defects, acting as a biological barrier preventing the invagination of soft tissue into the place where bone repair is taking place.

In implantodontics or periodontics, it can be used in the treatment of vertical and horizontal periodontal defects, lifting the maxillary sinus floor, in pre- and peri-implant surgery, in covering exposed implants or in covering bone defects to form new bone in alveolar edges. Pratix is not designed for any use other than that indicated.

#### CONTENTS AND PRESENTATION

Membrana Biológica Sintética – Absorvível (Pratix) is presented STERILE in rectangular and disc forms, packaged in a unitary form. The membrane is sold in: aluminum packaging (primary medical grade packaging) thermally sealed, packed in cardboard box - final packaging. The box is labeled with the relevant data and information.

The instructions for using the medical product (instructions for use) are available in the final packaging, one of which must be fixed in the protocol and the other delivered to the patient.

It is the responsibility of the hospital institution to ensure the identification and traceability of the medical product at the time of use through the adhesive label provided.

Figure 2 - Detail of the primary packaging in the final form in which Pratix is marketed. The yellow chemical indicator pack (with the symbol in the center) indicates that it has not yet been sterilized. The packaging containing the red indicator indicates that the material has already been irradiated.

The product codes Membrana Biológica Sintética – Absorvível, model Pratix in its rectangular shape are presented in the following tables:

MEMBRANA BIOLÓGICA SINTÉTICA - ABSORVÍVEL Pratix (Rectangular) 1 unit per package		
Code	Size	Thickness
992.R.150.1520	15 x 20mm	150 - 200 μ
992.R.150.1525	15 x 25mm	150 - 200 μ
992.R.150.2025	20 x 25mm	150 - 200 μ
992.R.150.3030	30 x 30mm	150 - 200 μ
992.R.250.1525	15 x 25mm	250 - 350 μ

The instructions for use and self-adhesive labels for product traceability are available inside the packaging.

#### SUPPORT MATERIAL

Supporting materials that come with the medical device are:

Instruction for using the product and identification labels for protocol and patient.

#### TRACEABILITY

Guidance to health professionals regarding the use and handling of traceability labels.

Self-adhesive labels (five) for product traceability are available inside the final packaging. The following traceability items are found on the labels: product reference code, manufacturing batch, manufacturer's logo and address, date of manufacture, validity, product registration number, and other information.

It is mandatory for the dentist or doctor to attach a label (label No. 1) to the patient's medical record or equivalent document, in order to comply with the current health legislation.

In order to continue the traceability cycle, the manufacturer recommends that the health professional assign the other available labels (labels n°2 to n°5) as follows:

Label Nº2: deliver to the patient

Label Nº3: fix in the tax documentation that generates collection (AIH in the case of a SUS patient, or invoice in the case of a patient assigned by the complementary health system).

Label Nº4: forward to supplier control

Label Nº5: forward to the control of the responsible (main) surgeon.

It is the sole responsibility of the dental hospital institution to ensure the identification and traceability of the medical product at the time of use, through the adhesive labels provided.

#### COMPOSITION

The product is composed form of poly (Lactic Acid) C6H8O4 and poly (Glycolic Acid) C4H4O4 in the proportion of 80/20%.

The chemical composition of the product and the semicrystalline structure were determined in the physico-chemical characterization studies and their results follow the standards of the American Society for Testing and Materials (ASTM) F1925 - 09 Standard Specification for Semi-Crystalline Poly (lactide) Polymer and Copolymer Resins for Surgical Implants, 1999.

In monitoring the physical and chemical characteristics of Pratix, the determined components are presented in the following percentages: pH = 6-8. Heavy metals: As = <10ppm ; Sb = <10ppm ; Cd = <10ppm ; Bi = <10ppm ; Hg = <10ppm ; Pb = <30ppm ; Sn = <10ppm and total heavy metals in the form of lead: <10ppm.

#### RECOMMENDATIONS

The patient must undergo a preoperative evaluation so that the best set of materials is determined, and the most appropriate techniques are established to improve the use of the product.

Only qualified professionals with knowledge of the surgical techniques and procedures necessary for the proper application and use of the product should use the product. Once the product is marketed in sterile form, the professional should be responsible for handling it following aseptic criteria.

The surgeon must ensure that the aseptic conditions are prepared and maintained in the surgical area and that the entire team is properly dressed. USE ONLY SPECIFIC AND STERILE SURGICAL INSTRUMENTS. Apply the product to viable and well-vascularized tissue.

The membrane must be cut beyond the defect margins to provide adequate protection to the defect and increase stability. However, be sure not to cut the membrane much larger than the defect, as this can cause interference with the flap's blood supply.

Important: Avoiding membrane exposure is the most frequent postoperative complication of GTR. If the membrane is exposed, strict plaque control is essential to prevent infection. Prepare thick flap. The blood supply of the flap depends on its thickness. There is a relationship between the thickness of the flap that covers the membrane and the occurrence of gingival recession after surgery.

For best results, leave approximately 1.5 mm of flap thickness to cover the membrane.

Ensure coronal positioning of the membrane. The more coronal the membrane is positioned, the more space for regeneration will be generated below the membrane. It is recommended to support the shape of the membrane by filing the area with bone grafts.

Smoking is an aggravating factor in periodontal disease; however, it is not an isolated factor, since smokers have poor plaque control.

#### USE INFORMATION

The use of the product with unsuitable surgical techniques and / or procedures and conditions may harm the patient, leading to unsatisfactory results.

Ideally, the surgeon should establish a post-operative control protocol including the measures taken in the situations encountered.

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It is also recommended that the surgeon keep a follow-up report of the evolution of the case together with the results. In the presence of tumors or infection the use of PratiX should be avoided.

The adhesive label accompanying the product must be attached to the patient's protocol to be kept in the clinic or hospital. Adopt classic prophylactic treatment described in the literature for immediate and late pre- and post-operative.

Like all biomaterial, it must be kept in places with low humidity and cool temperature. Technical and operational preparation for product use:

- The surgeon should check the integrity of the packaging - Do not use if the packaging is broken.
- The surgeon should check the expiration date of the product - Do not use if it is out of date.
- The surgical team should be responsible for preparing the patient and the area to perform a sterile procedure.
- The surgical team should be responsible for preparing sterile instruments and support materials (including saline to hydrate products).
- The team should be responsible for observing the patient's and team's attire for sterile procedure.

#### INSTRUCTIONS FOR USE

1. Carefully remove the membrane from the packaging using sterile gloves and instruments.
2. Cut the membrane so that it fits the bone defect in the best possible size to protect the defect, considering a maximum of 3 to 5 mm beyond the edges of the bone defect.
3. Handle and keep the product in a sterile field throughout the procedure (from opening, adjusting the size of the membrane, placing it in the surgical site until closing and finalization of the procedure).
4. Choose the most suitable membrane for the defect area in terms of size and shape at 5 mm beyond the defect margin. **Attention:** excess membrane interferes with the blood supply to the flap. The margin of the membrane must adapt close to the surrounding healthy bone tissue.
5. In most cases, the size of the membrane is adjusted by cutting out the lateral and apical portion. The collar portion should be left intact as much as possible.

#### POST-OPERATIVE CARE

During the healing period, the patient should use 0.1% chlorhexidine solution for mouthwash twice a day and apply the gel at 1%. Instruct the patient to start tooth brushing with extra soft bristle brushes 1 week after surgery. Instruct the patient to return weekly for dental cleaning until complete healing of the surgical wound (usually 4 to 6 weeks).

**Attention:** mechanical cleaning should not be restarted early in the defect area. In alcoholic patients (alcoholics), these membranes can be reabsorbed in less than 8 to 10 weeks for normal patients, since alcohol acts as an organic solvent of PLGA (lactic acid-co-glycolic acid) PLGA.

#### CONTRAINDICATIONS

Due to the lack of specific studies, the use of this product is contraindicated in patients with the following clinical signs / symptoms:

- Uncontrolled diabetes;
- Hypercalcemia
- General renal dysfunction;
- Severe vascular or neurological disease
- Severe liver disease;
- Patients on high-dose corticosteroid therapy;
- vascular damage in the implant area;
- Osteomyelitis in the area of surgery;
- Patients with generalized infection;
- Pregnant and lactating women (at medical discretion);
- Non-cooperative patients who are unable or unable to follow postoperative instructions, including those who abuse drugs and / or alcohol.

In the presence of tumors, the use of the product associated with biomaterials of an organic nature is contraindicated. Smoking is a major aggravating factor in periodontal disease, however, it is not an isolated factor, but smokers have poor plaque control.

In alcoholic patients (alcoholics), these membranes can be reabsorbed in less than 8 to 10 weeks for normal patients, since alcohol acts as an organic solvent of PLGA (lactic acid-co-glycolic acid) PLGA.

#### POSSIBLE ADVERSE EFFECTS

Possible complications such as swelling, infections, pain and complications associated with anesthesia may occur because they are common to implant surgery. As with any type of implant surgery, the patient may experience discomfort in the first postoperative days.

In addition to the aforementioned, normal reactions in periodontal surgeries may be found such as: swelling of intraoral tissue, thermal sensitivity, gum recession, some loss of height of crystal bone, etc. Spontaneous exfoliation may occur if Synthetic Membrana Biológica Sintética - Absorbível has not been properly covered with gingival tissue. In the clinical studies evaluating the Synthetic Biological Membrane - Absorbible, no immune reactions or hypersensitivity reactions were found.

#### SPECIAL PRECAUTIONS AND CLARIFICATIONS ABOUT STORAGE AND TRANSPORTATION OF THE PRODUCT

The packaging must be intact upon receipt. Do not use the medical product if the packaging is broken. The product must be transported and stored at room temperature. Avoid humidity and excessive heat.

The storage location of the medical device must be clean, dry and bright in order to maintain the ideal conditions for storage and transport, as well as its physical and chemical integrity.

The carrier must be informed of the content and delivery time. The precautions to be taken during transport are shown on the packaging (as shown in table 1), thus ensuring protection of the medical product from dispatch to delivery to the customer. Do not use the product after the expiration date.

Keep the product away from children.

SÍMBOLO	MEANING
	Single-use product.
	Caution, consult accompanying documents.
	Fragile, handle with care.
	Maximum temperature limit.
	Keep dry.
	Date of fabrication.
	Manufacturer
	Valid until
	Lot code.
	Catalog number.
	Sterilized by Gamma Ray.
	Do not use if packaging is damaged.

Table 1 - Information related to use, safety, handling, transport, storage and traceability of the product Synthetic Biological Membrane - Absorbable (PratiX) were placed in the form of symbols.

ISO 15223-1: 2012 Medical devices - symbols to be used with medical device labels, labeling and information to be supplied - part 1: general requirements; ISO 15223-2: 2010 / Medical devices - symbols to be used with medical device labels, labeling, and information to be supplied - part 2: symbol development, selection and validation.

#### INSTRUCTIONS AND PROCEDURES FOR DISPOSAL

The medical products described in this instruction for use are for SINGLE USE, so what is not used (leftovers) must be appropriately nich packaged and discarded by the hospital or dental institution. The leftovers should be autoclaved before the final de-characterization by applying a disposal cycle or with equivalent lethality (eg 121 ° C for 30 minutes or 132-134 ° C 12-10 minutes).

The product whose packaging is damaged or tampered with (compromising the sterility of the product) or whose expiration date must NOT be used. The same procedures used to dispose of product leftovers must be followed. BAUMER recommends that products in the situations described above be mechanically deformed, by making the remains unusable by cutting or tearing, after which they must be clearly and visibly identified, indicating their situation, that is, UNFIT TO USE. The methods and procedures for disposing of the medical product used must ensure its complete lack of character, preventing any possibility of reuse.

The mischaracterization of the medical product is the sole responsibility of the medical surgeon or dentist using the product as well as the disposal methods and procedures used.

#### PRECAUTIONS AND WARNINGS

The use of the product without surgical techniques and/or procedures and/or appropriate asepsis conditions, can harm the patient leading to unsatisfactory results.

In order not to expose the membrane and avoid a post-operative inflammatory response, approximately 1.5 mm of flap thickness is required to cover the membrane. In addition, if there is tension in the place where the flap covers the membrane, blood supply may be impeded in the case of very thin flaps. The coronal position of the membrane is an important factor in predicting the amount of tissue regeneration, the more coronal the membrane is positioned, the greater the regeneration achieved. Studies have confirmed that when the membrane collapsed within the area of the intraosseous defect, less regeneration occurred, so if the membrane is placed more coronal, more space for regeneration will be acquired, however, it is very useful to support the membrane shape with grafts bone.

One cause of early exposure of the membrane is necrosis of the flap that covers it. In guided bone regeneration, the blood supply of the flap depends on its thickness, since the supply that comes from the bone is impeded by the membrane. Therefore, the gingival thickness of the flap is crucial. There is a relationship between the thickness of the flap that covers the membrane and the occurrence of gingival recession after surgery.

Smokers, drug users, alcoholics or malnourished patients should be warned about the increased incidence of bone non-union;

#### STABILITY

Material for single use. Not re-sterilizable. Prohibited reprocessing. Do not use if the packaging is broken. Shelf life indicated on the box.

#### STERILITY

Product marketed STERILE. Sterilized with GAMMA RADIATION. NON-STERILIZABLE MATERIAL.

#### WARRANTY

The guarantee will be applied to medical products as long as it is used in accordance with the previous sections of this instruction for use: FUNDAMENTS OF PRODUCT ACTION, CONTENTS AND PRESENTATION, SUPPORT MATERIALS, TRACEABILITY, COMPOSITION, INDICATION, USE INFORMATION, INSTRUCTIONS FOR USE, POST OPERATIVE CARE, CONTRAINDICATION, POSSIBLE ADVERSE EFFECTS, SPECIAL PRECAUTIONS AND CLARIFICATIONS ABOUT STORAGE AND TRANSPORTATION OF THE PRODUCT , INSTRUCTIONS AND PROCEDURE FOR DISPOSAL, PRECAUTIONS AND WARNINGS, STABILITY AND STERILITY.

#### COMPLAINT

Any customer or user of these medical products who wishes to claim, if he is not satisfied with the services and / or products offered by BAUMER, such as: identification, traceability, security, effectiveness, performance, should contact the authorized BAUMER distributor.

Any defect or suspected defect, or problems with related medical products, must be reported immediately to BAUMER or your authorized BAUMER distributor.

Any serious adverse effects that affect the health or safety of the patient or user, such as malfunction or damage to medical products, serious problems or death caused by, or associated with, the use of these products, must be reported to BAUMER.

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