# Genius

Techgraft

**Membranas Regeneradoras**

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MEMBRANA BIOLÓGICA ACELULAR DE ORIGEM BOVINA		
TECHGRAFT		
1 unidade por embalagem		
Código	Tamanho	Espessura
985.R.1010	1,0 x 1,0 cm	150-250m
985.R.1015	1,0 x 1,5 cm	150-250m
985.R.1020	1,0 x 2,0 cm	150-250m
985.R.1025	1,0 x 2,5 cm	150-250m
985.R.1030	1,0 x 3,0 cm	150-250m
985.R.1035	1,0 x 3,5 cm	150-250m
985.R.1040	1,0 x 4,0 cm	150-250m
985.R.1045	1,0 x 4,5 cm	150-250m
985.R.1050	1,0 x 5,0 cm	150-250m
985.R.1055	1,0 x 5,5 cm	150-250m
985.R.1060	1,0 x 6,0 cm	150-250m
985.R.1515	1,5 x 1,5 cm	150-250m
985.R.1520	1,5 x 2,0 cm	150-250m
985.R.1525	1,5 x 2,5 cm	150-250m
985.R.1530	1,5 x 3,0 cm	150-250m
985.R.1535	1,5 x 3,5 cm	150-250m
985.R.1540	1,5 x 4,0 cm	150-250m
985.R.1545	1,5 x 4,5 cm	150-250m
985.R.1550	1,5 x 5,0 cm	150-250m
985.R.1555	1,5 x 5,5 cm	150-250m
985.R.1560	1,5 x 6,0 cm	150-250m
985.R.2020	2,0 x 2,0 cm	150-250m
985.R.2025	2,0 x 2,5 cm	150-250m
985.R.2030	2,0 x 3,0 cm	150-250m
985.R.2035	2,0 x 3,5 cm	150-250m
985.R.2040	2,0 x 4,0 cm	150-250m
985.R.2045	2,0 x 4,5 cm	150-250m
985.R.2050	2,0 x 5,0 cm	150-250m
985.R.2055	2,0 x 5,5 cm	150-250m
985.R.2060	2,0 x 6,0 cm	150-250m
985.R.3030	3,0 x 3,0 cm	150-250m
985.R.3035	3,0 x 3,5 cm	150-250m
985.R.3040	3,0 x 4,0 cm	150-250m
985.R.3045	3,0 x 4,5 cm	150-250m
985.R.3050	3,0 x 5,0 cm	150-250m
985.R.3055	3,0 x 5,5 cm	150-250m
985.R.3060	3,0 x 6,0 cm	150-250m
985.R.4040	4,0 x 4,0 cm	150-250m
985.R.4045	4,0 x 4,5 cm	150-250m
985.R.4050	4,0 x 5,0 cm	150-250m
985.R.4055	4,0 x 5,5 cm	150-250m
985.R.4060	4,0 x 6,0 cm	150-250m
985.R.5050	5,0 x 5,0 cm	150-250m
985.R.5055	5,0 x 5,5 cm	150-250m
985.R.5060	5,0 x 6,0 cm	150-250m
985.R.6060	6,0 x 6,0 cm	150-250m

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**Techgraft**

**DESCRIPTION**

TECHGRAFT is a biological acelular membrane of bovine origin, which is implantable and resorbable, it is designed to be used as a biological barrier for tissue regeneration surgical procedure in the field of dentistry.

TECHGRAFT is made from bovine pericardium, and the end product is a natural high purity collagen membrane which acts as a biological barrier when placed between the grafting material and the soft tissue, enabling the guided regeneration of the bone and connective tissues as well.

In the manufacture of TECHGRAFT, the tissue of origin subjected to a standardized and controlled production. One of the particularities of this multiple process is the preservation of the characteristics of the type I collagen present in the tissue of origin. The end product is biocompatible, with high purity, free of heavy and biological contaminant, acelular, noncytotoxic, nonimmunogenic, nonpyrogenic. Its porosity favors the cell repopulation and the neovascularization, essential conditions for the remodeling process and for the reconstruction of the bone and conjunctive tissues.

TECHGRAFT is offered in a rectangular membrane shape or in different sizes of discs.

TECHGRAFT is provided STERILE. It is sterilized with GAMMA RAYS.

TECHGRAFT is a SINGLE USE PRODUCT, therefore, REPROCESSING IS PROHIBITED.

**PERFORMANCE AND ACTION FOUNDATIONS**

The usage of biological barriers as TECHGRAFT to repair the soft tissue represents an alternative to autogenous grafting of the soft tissues. The product is ready to use and it eliminates the necessity of removing tissues from the patient, reducing the surgical time and the possibility of complications to happen, minimizing pains and morbidities associated to the procedure. The high purity of TECHGRAFT allows a naturaland harmonious incorporation to the patient’s tissue, creating the ideal coverage of the recession or the effective regeneration of the keratinous tissue.

TECHGRAFT works as a fiscal battery on Guided Tissue Regeneration (GTR), associated or not with synthetic or natural bone.

The concept is the usage of a passive barrier to prevent the unwelcome migration of epithelial tissue and gingival conjunctive tissue to the bone defect during the remodeling.

This membrane’s flexibility is fundamental to its adaption to the anatomic structures where they are going to be placed during the surgery.

The absorption speed of the membrane may be influenced by local characteristics, as vascularisation of the area, thickness of the mucos, inflammation level consequential of the surgical act, modifications on the pH of the tissue, local infections and partial exposures to the mouth. The results obtained in studies with the product and literature indicate partial absorption in 30 days, and complete in 60 days.

**PRODUCT COMPOSITION**

The analysis of the qualitative and quantitative composition of phase indicated that Biological Acelular Bovine-derived Membrane (TECHGRAFT) is a biological material composed of collagen, especially of the type I.

The chemical composition of the product and the semi crystalline structure were determined in physico-chemical characterization studies using scanning calorimetry (DSC) and mass spectrometry with inductively coupled plasma (ICP). The tests and acceptance criteria follow the standards of the American Society for Testing and Materials (ASTM) F2212-11 Standard Guide for Characterization of Type I Collagen as Starting Material for Surgical Implants and Substrates for Tissue Engineered Medical Products (TEMPS), 2011.

On monitoring of physico-chemical characteristics of the product, the components are presented in the following percentages: pH = 6-8. Heavy metals: As <= <10 ppm ; Sb <= <10 ppm; Cd <= <10 ppm; Bi = <10 ppm ; Hg <= <10 ppm ; Pb = < 30 ppm; Sn = < 10 ppm and total heavy metals in the form of lead: <10ppm Collagen concentration: 13.1 g/100.

Concentration of amino acids: cysteine 25.1 g/100, hydroxyproline 13.7 g/100 tryptophan 16.3 g/100 and tyrosine 15.8 g/100.

**PRESENTATIONS OF THE MEDICAL DEVICE**

Biological Acelular Membrane of Bovine Origin (TECHGRAFT) is presented STERILE in square, rectangular and circular forms.

The mem-brane is marketed in: blister pack of virgin polyethylene terephthalate (PET) medical degree and thermally sealed aluminium envelope packaged in cardboard box along with the instructions for use and the traceability labels. The box is sealed and labeled with the data and information pertinent to the product.

The instruction for use of the medical device is provided within the final packing-box.

The total of five self-adhesive labels with product information are available in the final packaging.

It is the responsibility of the hospital or dental institution to ensure the identification and traceability of the medical product at the time of use through the adhesive labels provided.

The TECHGRAF product codes in their rectangular form are shown in the following table:

BIOLOGICAL ACELLULAR MEMBRANE OF BOVINE ORIGIN		
TECHGRAFT		
1 unit per carton		
Code	Diameter	Thickness
985.R.1010	1,0 x 1,0 cm	150-250m
985.R.1015	1,0 x 1,5 cm	150-250m
985.R.1020	1,0 x 2,0 cm	150-250m
985.R.1025	1,0 x 2,5 cm	150-250m
985.R.1030	1,0 x 3,0 cm	150-250m
985.R.1035	1,0 x 3,5 cm	150-250m
985.R.1040	1,0 x 4,0 cm	150-250m
985.R.1045	1,0 x 4,5 cm	150-250m
985.R.1050	1,0 x 5,0 cm	150-250m
985.R.1055	1,0 x 5,5 cm	150-250m
985.R.1060	1,0 x 6,0 cm	150-250m
985.R.1515	1,5 x 1,5 cm	150-250m
985.R.1520	1,5 x 2,0 cm	150-250m
985.R.1525	1,5 x 2,5 cm	150-250m
985.R.1530	1,5 x 3,0 cm	150-250m
985.R.1535	1,5 x 3,5 cm	150-250m
985.R.1540	1,5 x 4,0 cm	150-250m
985.R.1545	1,5 x 4,5 cm	150-250m
985.R.1550	1,5 x 5,0 cm	150-250m
985.R.1555	1,5 x 5,5 cm	150-250m
985.R.1560	1,5 x 6,0 cm	150-250m
985.R.2020	2,0 x 2,0 cm	150-250m
985.R.2025	2,0 x 2,5 cm	150-250m
985.R.2030	2,0 x 3,0 cm	150-250m
985.R.2035	2,0 x 3,5 cm	150-250m
985.R.2040	2,0 x 4,0 cm	150-250m
985.R.2045	2,0 x 4,5 cm	150-250m
985.R.2050	2,0 x 5,0 cm	150-250m
985.R.2055	2,0 x 5,5 cm	150-250m
985.R.2060	2,0 x 6,0 cm	150-250m
985.R.3030	3,0 x 3,0 cm	150-250m
985.R.3035	3,0 x 3,5 cm	150-250m
985.R.3040	3,0 x 4,0 cm	150-250m
985.R.3045	3,0 x 4,5 cm	150-250m
985.R.3050	3,0 x 5,0 cm	150-250m
985.R.3055	3,0 x 5,5 cm	150-250m
985.R.3060	3,0 x 6,0 cm	150-250m
985.R.4040	4,0 x 4,0 cm	150-250m
985.R.4045	4,0 x 4,5 cm	150-250m
985.R.4050	4,0 x 5,0 cm	150-250m
985.R.4055	4,0 x 5,5 cm	150-250m
985.R.4060	4,0 x 6,0 cm	150-250m
985.R.5050	5,0 x 5,0 cm	150-250m
985.R.5055	5,0 x 5,5 cm	150-250m
985.R.5060	5,0 x 6,0 cm	150-250m
985.R.6060	6,0 x 6,0 cm	150-250m

The TECHGRAF product codes in their circular form (discs) are presented in the following table:

BIOLOGICAL ACELLULAR MEMBRANE OF BOVINE ORIGIN		
TECHGRAFT		
1 unit per carton		
Code	Diameter	Thickness
985.D.10	1,0 cm	150-250m
985.D.15	1,5 cm	150-250m
985.D.20	2,0 cm	150-250m
985.D.25	2,5 cm	150-250m
985.D.30	3,0 cm	150-250m
985.D.35	3,5 cm	150-250m
985.D.40	4,0 cm	150-250m
985.D.45	4,5 cm	150-250m
985.D.50	5,0 cm	150-250m
985.D.55	5,5 cm	150-250m
985.D.60	6,0 cm	150-250m
985.D.65	6,0 cm	150-250m

**SUPPORT MATERIAL**

The support materials that accompany the medical product are: Product instruction for use and traceability labels.

**TRACEABILITY –** guidelines to the health professional regarding the use and handling of the traceability labels.

Self-adhesive labels (five) for product traceability are available inside the final packaging. The following information and traceability items are found on the labels:

Model, technical name and trade name of the product

Logo and address of the manufacturer,

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## INDICATIONS

TECHGRAFT is indicated to be used as an absorbable biological barrier with physical function of separation and protection against migration or unwanted cell invasion in tissue regeneration dental surgeries.

TECHGRAFT can be used in association with Synthetic or natural bone grafts, including the treatment of periodontal diseases with or without bone loss.

The product is indicated to be used in dental surgeries as a resorbable material to: oral and maxillofacial reconstructions involving gingival retraction, increase alveolar ridge (bone and soft tissue), reconstruction of alveolar ridge for prosthetic treatment, increase around implants with late or immediate placement after extraction, targeted ridge augmentation to late implant, filling of bone defects after recession of roots, removal of cysts or retained teeth surgery, guided regeneration procedures in periodontal defects and the preservation of the extraction place.

Biological Acellular Membrane of Bovine Origin (TECHGRAFT) is not designed for other uses other than those indicated.

## USAGE DETAILS

Only qualified professionals and with knowledge of surgical techniques and the needed procedures for the choice, application and appropriate use of the product should make use of it.

The use of the product with inappropriate surgical techniques and conditions may impair the patient, leading to unsatisfactory results.

The patient should undergo a pre-operative evaluation so that the best set of materials is determined and the appropriate techniques are established for the use of the product.

Ideally, the surgeon should establish a post-operative control protocol including the measures taken in the situations encountered. It is recommended that the surgeon also hold a follow-up report on the evolution of the case along with the results.

The surgeon should provide the classic prophylactic treatment described in the literature for immediate and late, pre and postoperative period.

The product should only be used for the purpose for which it is intended.

The membrane may be cut for better adaptation to the defect. This step of TECHGRAFT adaptation to the shape of the defect can be done before or after the hydration in saline.

Before the implantation, the membrane should be hydrated in sterile saline solution.

After this stage, the membrane beco-mes more transparent, but still keeping its white color, and flexible. Hydration helps on the handling and does not interfere in the physical chemical nor structural characteristics of the membrane.

## INSTRUCTIONS FOR USE:

a – At room temperature, put approximately 50 ml of sterile saline solution in a sterile container.

b – Check if the packaging is intact. (Make sure that the package has not been tampered with and that the product is within the shelf life)

c – Remove the product from the primary packaging with a sterile clamp, preferably with flattened pit and put it in the container with saline.

d –With the assistance of tweezers, keep the product immersed in the saline solution. Once the product remain submerged without the pressure of the clamp, leave it in the solution until it reaches a flexible and homogeneous transparent aspect (from 5 to 20 minutes).

e – Cut the product into the desired shape and adapt on the imperfection. TECHGRAFT has a good performance in relation to modeling.

## CONTRAINDICATIONS

Due to the lack of specific studies, the use of this product is contraindicated in patients with the following clinical conditions:

- Uncontrolled Diabetes;
- General renal dysfunction;
- Vascular disease or severe neurological disorder;
- Severe liver diseases;
- Patients under high-dose corticotherapy;
- Vascular damage in the area of the implant;
- Osteomyelitis in the surgery area;
- Patients with sepsis;
- Pregnant women and lactating mothers (in the discretion of the doctor);
- Uncooperative patients who will not or cannot follow the postoperative instructions, including those who abusively use drugs and/or drink alcohol.

In case there are tumors, the use of the product associated with the organic nature biomaterials is contraindicated. Smoking is a major aggravating factor of the periodontal disease, however, is not an isolated factor. Smokers have low dental plaque control.

## POSSIBLE ADVERSE EFFECTS

Possible complications as swelling, infection, pain and complications associated with anesthesia may occur for being common in implant surgeries. As any type of implant surgery, the patient may experience discomfort in the early days of the postoperative period.

In addition to the aforementioned, normal reactions in periodontal surgery may be found as: intraoral tissue swelling, heat sensitivity, gum recession, some loss of bone height, etc. Spontaneous exfoliation may occur if the product has not been properly covered with gum tissue.

In clinical studies of product evaluation, it has not been found immune reactions or hypersensitivity reactions.

## SPECIAL PRECAUTIONS AND CLARIFICATIONS ABOUT THE STORAGE AND TRANSPORTATION OF THE MEDICAL DEVICE

Medical products described in this technical report are for SINGLE USE, therefore what is not used (leftovers) should be adequately destroyed and discarded by the hospital or dental institution. Autoclaving should be done with the leftover before the final distortion by making a disposal cycle or something with equivalent lethality (e.g. 121° C for 30 minutes or 132-134° C for 12-10 minutes).

The product which the packaging has been damaged or tampered with (compromising the sterility of the product) or with expired period of validity **SHOULD NOT BE USED**. The same procedures employed to dispose the leftover product should be followed. BAUMER recommends the products in the situations described above to be mechanically deformed through the destruction of the remains by cutting or ripping, and then they should be identified clearly and visibly, indicating its situation, in other words, **INAPPROPRIATE FOR USE**.

The methods and procedures for the disposal of the medical product used should ensure the complete disfigurement of it, preventing any possibility of its re-use. The responsibility for the disfigurement of the medical device is solely on the surgeon or dentist who are using the product as well as for the disposal methods and procedures used.

SYMBOL	MEANING
	Single use product.
	Caution, consult accompanying documents.
	Fragile, handle carefully.
	Upper temperature limit.
	Keep dry.
	Manufacturing date.
	Manufacturer
	Expiring date.
LOT	Batch Code.
REF	Catalogue number.
	Sterilized with gamma rays.
	Do not use if the package is damaged.

Tabela 1 - The information concerning use, safety, handling, transportation, storage and traceability of the product Bovine Pericardium Membrane were placed in the form of symbols.

**Source:** ISO 15223-1:2012 Medical devices-symbols to be used with medical device labels, labelling and information to be supplied-part 1: general requirements; ISO 15223-2:2010/Medical devices-symbols to be used with medical device labels, labelling and information to be supplied- part 2: symbol development, selection and validation.

## INSTRUCTIONS AND PROCESURES OF DISPOSAL OF THE MEDICAL DEVICE

Medical products described in this technical report are for SINGLE USE, therefore what is not used (leftovers) should be adequately destroyed and discarded by the hospital or dental institution. Autoclaving should be done with the leftover before the final distortion by making a disposal cycle or something with equivalent lethality (e.g. 121° C for 30 minutes or 132-134° C for 12-10 minutes).

The product which the packaging has been damaged or tampered with (compromising the sterility of the product) or with expired period of validity **SHOULD NOT BE USED**. The same procedures employed to dispose the leftover product should be followed.

BAUMER recommends the products in the situations described above to be mechanically deformed through the destruction of the remains by cutting or ripping, and then they should be identified clearly and visibly, indicating its situation, in other words, **INAPPROPRIATE FOR USE**.

The methods and procedures for the disposal of the medical product used should ensure the complete disfigurement of it, preventing any possibility of its re-use. The responsibility for the disfigurement of the medical device is solely on the surgeon or dentist who are using the product as well as for the disposal methods and procedures used.

## WARNINGS AND PRECAUTIONS

The use of the product without surgical techniques and/or procedures and/or adequate asepsis conditions can harm the patient and lead to unsatisfactory results.

Not to accur membrane exposure and prevent postoperative inflammatory response, it is necessary approximately 1.5 mm of thickness of the flap to cover the membrane. In addition, if there is tension in the place where the flap covers the membrane, the blood supply can be prevented in case the flap isf very slender.

The surgeon should completely know the medical and surgical aspects of the medical device.

Smoker, drug users, alcoholics or malnourished patients should be alerted about the increased incidence of nonunion bone.

## WARRANTY

The warranty will be applied to the provided medicinal products that are used according to the established in the previous sections of this instruction for use: SUPPORT MATERIALS, INDICATIONS, USAGE INFORMATION, CONTRAINDICATIONS, CARE AND PRECAUTIONS, WARNINGS, SPECIAL PRECAUTIONS AND CLARIFICATIONS ON THE STORAGE AND TRANSPORT OF THE MEDICAL DEVICE AND STERILITY.

## COMPLAINT

All client or user of these medical products who would like to complain, if not satisfied with the services and/or products offered by BAUMER, such as: Identification, traceability, safety, efficiency, performance, should contact the authorized BAUMER distributor. Any defect or suspicion of defect, or problems related to the medical products must be immediately reported to BAUMER or to the authorized BAUMER distributor.

In the occurrence of any adverse effect related to the product that affects the health or safety of the patient or user, the professional dentist or doctor should notify the competent sanitary authority (ANVISA) through the website http://portal.anvisa.gov.br and report what happened to BAUMER or the authorized BAUMER Distributor.

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