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#### MATERIAL SUPPORT

The support materials that accompany the medical product are:  
- Instructions for use of the product  
- Identification tag for protocol and patient.

#### INDICATIONS

The support Enxerto Ósseo Cerámico Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) is indicated as a bone substitute with osteoconductive function presented in the form of granules or blocks indicated as filling material for surgical bone cavities where it is needed later replacement by bone tissue. The product can be used alone or combined with other kind of autologous bone grafts, autogenous or animal origin or synthetic grafting procedures in the various areas of Periodontology, Implantology, Surgery Periapical surgery and maxillo-facial, and can also be applied in areas of orthopedics, traumatology, neurosurgery, plastic surgery or repair, cranial surgery and oncology.

#### USE OF INFORMATION

Only qualified professionals with knowledge and surgical techniques and procedures for selection, application and appropriate use of the product must make use of it.  
The use of the product with non-surgical techniques and appropriate conditions may harm the patient leading to unsatisfactory results. The patient should undergo preoperative evaluation for the best set of material is determined and the most appropriate techniques are established for use of the product.  
Ideally the surgeon should establish a protocol for control of postoperative including measures taken in the situations encountered. It is also recommended that the surgeon maintains a monitoring report on progress of the case together with the results. The surgeon should take prophylactic treatment classically described in the literature for pre-and postoperative early and late. The adhesive label with the product should be fixed in the protocol the patient being kept in the clinic or hospital.

#### TECHNICAL AND OPERATIONAL PREPARATION FOR PRODUCT USE

The patient should undergo preoperative evaluation for the best set of material is determined and the most appropriate techniques are established for use of the product.

The surgeon should take prophylactic treatment classically described in the literature for pre-and postoperative early and late. The surgeon should check the integrity of the package - Do not use if pouch is damaged.  
The surgeon should check the validity of the product - Do not use if you are out of date.  
The surgical team will be responsible for preparing the patient and the area to perform sterile procedure.  
The surgical team will be responsible for preparing instruments and materials to support sterile (including saline for hydration products).  
The team will be responsible for observing the patient and the scrub team to sterile procedure.  
Check the integrity and validity of the product before opening.  
Using sterile technique.

#### INSTRUCTIONS FOR USE

1. Before opening the package should be prepared container for hydration of the Enxerto Ósseo Cerámico Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) and combination with other graft materials. The container should be sterile and fluid used for hydration. It is suggested to use sterile saline or the patient's blood taken under sterile conditions at the time of surgery.  
2. Observe the site preparation and proper use of sterile scrub of the physician and patient adviser to sterilize procedure.  
3. Check the integrity and validity of the product before opening.  
4. Using sterile technique. Remove the bottle and secondary end of the package.  
5. Remove the seal and the rubber bung.  
Warning: Carefully remove the seal, the seal can stick metal gloves.  
6. Add the contents of the vial to a sterile bowl and mix with sterile saline.  
7. The product Enxerto Ósseo Cerámico Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) can be used alone or in combination with other graft materials for bone remodeling such as autograft, freeze-dried demineralized bone matrix (DFDB human or xenogenic), deproteinized bovine bone, etc. Mix sterile saline solution to the graft with a sterile spatula.  
8. Warning: Before applying the product debríandona prepare the area of fibrous tissue, necrotic or fibrous and in the case of direct contact with bone decortication perform in order to provide direct contact with bone graft viable.  
9. Irrigate the surgical site before placing the graft. For best results, do any manipulation of the site before placing the graft.  
10. Use the product immediately after reconstitution with saline. With a sterile spatula to fill the bone defect.  
11. **Warning:** make sure there is adequate for hemostasis assuregarse que the material is retained in the graft site.  
12. Proceed as usual to close or cover the area deployed.

**Figura 2** - Details of the product Enxerto Ósseo Cerámico Bifásico (Models: GenPhos HA TCP

OrthoCer HA TCP, GenPhos XP and OrthoCer XP) being hydrated in saline immediately prior to application to the bone defect. The hydrated product should be mixed with a spatula and applied in failure.

#### INSTRUCTIONS AND PROCEDURES FOR DISPOSAL OF THE MEDICAL PRODUCT

The product Enxerto Ósseo Cerámico Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) single use, so that it not used (surplus) must be properly edescartado mischaracterized by the hospital or dental user. Should be carried out prior to autoclaving the remains of the disfigurement by applying a final cycle of disposal or another cycle with equivalent lethality (eg 121 ° C for 30 minutes or 12-10 minutes 123-134 ° C).  
The product whose packaging is damaged or inviolate (compromising the sterility of the product) or with expired shelf life should not be used. Should be followed the same procedures used for disposal of leftover product.  
The products BAUMER recommends that the situations described above are mechanically deformed by the destruction of the remains by the destruction of the glass bottle, cut or tear the packaging primary, secondary and final (box) of which must then be clearly identified and visible, indicating their status, ie, IMPROPER USE.  
The methods and procedures for disposal of medical products used must ensure the complete characterization of it, preventing any possibility of reuse. The mischaracterization of the medical product is the sole responsibility of the surgeon doctor or dentist that you are using the product as well as the methods and procedures used for disposal.

#### PRECAUTIONS AND CAUTIONS

For safety and effectiveness in use of product family Enxerto Ósseo Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) the following precautions and care should be observed:  
Only qualified professionals with knowledge and surgical techniques and procedures for proper selection and application should make use of the product.  
The surgeon doctor or dentist should be familiar with the procedures and recommendations for the use of GenPhos associated with other graft materials (see INFORMATION FOR USE).

Should only be used specific surgical instruments sterilized.  
The use of surgical instruments under conditions of preparation and sterilization of instruments and area, and vestments of the team.  
As with all surgical procedures should be taken care in the treatment of medically compromised patients and patients on steroid therapy for a long time, or patients taking anticoagulants.  
The family of products Enxerto Ósseo Bifásico (Models: GenPhos HA TCP, OrthoCer HA TCP, GenPhos XP and OrthoCer XP) is composed of single-use products that should not be re-sterilized. Open products should be discarded.

#### WARNING

The surgeon must learn is full of medical and surgical aspects of medical product.  
Smokers, drug users, alcoholics or malnutrition, should be warned about the increased incidence of non-union bone.

#### POSSIBLE ADVERSE EFFECTS

There were no reported adverse reactions to the product GenPhos HA TCP.  
Possible complications such as swelling, infection, pain and complications associated with anesthesia can occur because they are common to implant surgeries. As with any type of implant surgery the patient may feel discomfort in the first days after surgery.

#### SPECIAL PRECAUTIONS AND CLARIFICATIONS ON THE STORAGE AND TRANSPORT MEDICAL PRODUCT

The packaging must be intact upon receipt.  
Do not open containers should be used (do not use the medical product if the packaging is broken);  
The product must be transported and stored at room temperature.  
Avoid moisture and excessive heat.  
The storage of medical products must be clean, dry and bright in order to maintain the ideal conditions for storage and transportation, as well as the identification, their physical and chemical;  
The carrier must be informed about the content and delivery date. Care should be taken in transport are shown on the packaging, thus ensuring protection of the medical device from the dispatch to delivery to the customer.

#### CONTRAINDICATIONS

The family of Enxerto Ósseo Cerámico Bifásico (GenPhos HA TCP and OrthoCer Ha TCP) was not designed for other use other than that indicated. Do not use the product in presence of any of the contraindications listed below.  
In patients with uncontrolled metabolic disorders or abnormal metabolism of calcium.  
In patients with malignant tumors.  
In patients with presence or recent history of systemic or localized infection (infected wounds in the area or near the implant)  
In patients with pictures of recent infection or immune disorder or untreated  
In patients with immunological abnormalities.  
In patients with immature bone or cranium incomplete.  
In patients with severe osteopenia or osteoporosis.  
In patients who received or are to receive chemo- or radiotherapy near / at the site of implantation.  
In areas where the surrounding bone is not viable or unable to support or anchor the graft.  
In patients with severe renal dysfunction or abnormalities.  
Due to lack of specific studies using Enxerto Ósseo Cerámico Bifásico (GenPhos HA TCP and OrthoCer Ha TCP) should be avoided in patients with the following tables:  
Renal widespread  
Severe liver disease  
Alcoholism or drug abuse;  
Pregnancy and lactancia (the physician's discretion);  
Fever or leucocitosenão is recommended to use TCP GenPhos HA biomateriales associated with organic nature.

#### STERILITY

The medical product Enxerto Ósseo Cerámico Bifásico is marketed sterile and should not be re-sterilized.  
Should not be re-used leftover product

#### WARRANTY

The guarantee will be applied to medical products since it is used in accordance with the provisions in the previous sections of this instruction for use: SUPPORT MATERIALS, INSTRUCTIONS, INFORMATION FOR USE, CONTRAINDICATIONS, WARNINGS AND PRECAUTIONS, WARNINGS, PRECAUTIONS AND CLARIFICATIONS ON storage and transportation of medical products and sterility.

#### CLAIM

Every client or user of medical products that you wish to complain, if not satisfied with the services and / or products offered by BAUMER, such as identification, traceability, safety, effectiveness, performance, should contact the authorized distributor BAUMER.  
Any defect or suspected defect, or problems related to medical products, must be immediately reported to the authorized distributor or BAUMER. Any serious adverse effects affecting the health or safety of the patient or user, such as damage or malfunction of medical devices, major problem or death caused by or associated with the use of these products must be reported to BAUMER.

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#### Baumer S.A.